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# Owner's Record

The model and serial numbers are located at the rear. Record these numbers in the space provided below. Refer to these numbers whenever you call upon your Sony dealer regarding this product.

Model No. \_\_\_\_\_

Serial No. \_\_\_\_\_

## WARNING

**To reduce the risk of fire or electric shock, do not expose this apparatus to rain or moisture.**

**To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only.**

**No modification of this equipment is allowed.**

**THIS APPARATUS MUST BE EARTHED.**

**This unit has no power switch. To disconnect the main power, unplug the power plug. When installing the unit, incorporate a readily accessible disconnect device in the fixed wiring, or connect the power plug to an easily accessible socket-outlet near the unit. Do not position the ME equipment where it is difficult to unplug the power plug. If a fault should occur during operation of the unit, operate the disconnect device to switch the power supply off, or unplug the power plug.**

### Symbols on the products

#### Consult the Instructions for Use

Follow the directions in the Instructions for Use for parts of the unit on which this mark appears.

This symbol indicates the manufacturer, and appears next to the manufacturer's name and address.

### For the customers in the U.S.A.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate this equipment.

All interface cables used to connect peripherals must be shielded in order to comply with the limits for a digital device pursuant to Subpart B of part 15 of FCC Rules.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

### For the customers in Canada

CAN ICES-3 (B)/NMB-3(B)


This unit has been certified according to Standard CAN/CSA-C22.2 No.60601-1.

### For the customers in the U.S.A. and Canada

Model AC-80MD is Non-Patient Equipment.

This unit can not be used in the vicinity of patients.

#### \* Patient Vicinity


### Important safeguards/notices for use in the medical environments

1. All the equipments connected to this unit shall be certified according to Standard IEC60601-1, IEC60950-1, IEC60065 or other IEC/ISO Standards applicable to the equipments.
2. Furthermore all configurations shall comply with the system standard IEC 60601-1-1. Everybody who connects additional equipment to the signal input part or signal output part configures a medical system, and is

therefore, responsible that the system complies with the requirements of the system standard IEC60601-1-1. If in doubt, consult the qualified service personnel.

This unit can not be used in the patient environment.

#### \* Patient Environment

3. Do not touch the patient simultaneously while you are contacting with this equipment.
4. For this particular equipment, all accessory equipment connected as noted above, must be connected to mains via an additional isolation transformer conforming with the construction requirements of IEC60601-1 and providing at least Basic Insulation.
5. The leakage current could increase when connected to other equipment.
6. This equipment generates, uses, and can radiate radio frequency energy. If it is not installed and used in accordance with the instruction manual, it may cause interference to other equipment. If this unit causes interference (which can be determined by unplugging the power cord from the unit), try these measures: Relocate the unit with respect to the susceptible equipment. Plug this unit and the

susceptible equipment into different branch circuit.

Consult your dealer. (According to standard IEC60601-1-2 and CISPR11, Class B, Group 1)

### **Note**

The accessibility of the plug part of the DC connector of this AC adaptor is taken into account in the description of the patient environment. This adaptor can be used in the patient environment because the plug part of the DC connector does not come in contact with human body while this adaptor is connected with Sony medical printers, recorders, etc.

**Important EMC notices for use in the medical environments**

- The AC-80MD needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this Instructions for Use.
- The portable and mobile RF communications equipment such as cellular phones can affect the AC-80MD.

**Warning**

The use of accessories and cables other than those specified, with the exception of replacement parts sold by Sony Corporation, may result in increased emissions or decreased immunity of the AC-80MD.

| Guidance and manufacturer’s declaration-electromagnetic emissions                                                                                                                    |            |                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The AC-80MD is intended for use in the electromagnetic environment specified below.<br>The customer or the user of the AC-80MD should assure that it is used in such an environment. |            |                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Emission test                                                                                                                                                                        | Compliance | Electromagnetic environment-guidance                                                                                                                                                                                                                                                                                                                                                                             |
| RF emissions<br><br>CISPR 11                                                                                                                                                         | Group 1    | The AC-80MD uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>The AC-80MD is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions<br><br>CISPR 11                                                                                                                                                         | Class B    |                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Harmonic emissions<br><br>IEC 61000-3-2                                                                                                                                              | Class D    |                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Voltage fluctuations/flicker emissions<br><br>IEC 61000-3-3                                                                                                                          | Complies   |                                                                                                                                                                                                                                                                                                                                                                                                                  |

**Warning**

If the AC-80MD should be used adjacent to or stacked with other equipment, it should be observed to verify normal operation in the configuration in which it will be used.

| Guidance and manufacturer's declaration - electromagnetic immunity                                                                                                                |                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The AC-80MD is intended for use in the electromagnetic environment specified below. The customer or the user of the AC-80MD should assure that it is used in such as environment. |                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                |
| Immunity test                                                                                                                                                                     | IEC 60601 test level                                                                                                                                                                                                                  | Compliance level                                                                                                                                                                                                                      | Electromagnetic environment-guidance                                                                                                                                                                                                                                           |
| Electrostatic discharge (ESD)<br><br>IEC 61000-4-2                                                                                                                                | $\pm 6$ kV contact<br><br>$\pm 8$ kV air                                                                                                                                                                                              | $\pm 6$ kV contact<br><br>$\pm 8$ kV air                                                                                                                                                                                              | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                  |
| Electrical fast transient/burst<br><br>IEC 61000-4-4                                                                                                                              | $\pm 2$ kV for power supply lines<br><br>$\pm 1$ kV for input/output lines                                                                                                                                                            | $\pm 2$ kV for power supply lines<br><br>$\pm 1$ kV for input/output lines                                                                                                                                                            | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                            |
| Surge<br><br>IEC 61000-4-5                                                                                                                                                        | $\pm 1$ kV line(s) to line(s)<br><br>$\pm 2$ kV line(s) to earth                                                                                                                                                                      | $\pm 1$ kV differential mode<br><br>$\pm 2$ kV common mode                                                                                                                                                                            | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                            |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br><br>IEC 61000-4-11                                                                        | $< 5\% U_T$<br>( $> 95\%$ dip in $U_T$ ) for 0.5 cycle<br><br>$40\% U_T$<br>( $60\%$ dip in $U_T$ ) for 5 cycles<br><br>$70\% U_T$<br>( $30\%$ dip in $U_T$ ) for 25 cycles<br><br>$< 5\% U_T$<br>( $> 95\%$ dip in $U_T$ ) for 5 sec | $< 5\% U_T$<br>( $> 95\%$ dip in $U_T$ ) for 0.5 cycle<br><br>$40\% U_T$<br>( $60\%$ dip in $U_T$ ) for 5 cycles<br><br>$70\% U_T$<br>( $30\%$ dip in $U_T$ ) for 25 cycles<br><br>$< 5\% U_T$<br>( $> 95\%$ dip in $U_T$ ) for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. If the user of the AC-80MD requires continued operation during power mains interruptions, it is recommended that the AC-80MD be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br><br>IEC 61000-4-8                                                                                                                    | 3 A/m                                                                                                                                                                                                                                 | 3 A/m                                                                                                                                                                                                                                 | Power frequency magnetic fields should be at least characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                       |
| NOTE: $U_T$ is the a.c. mains voltage prior to application of the test level.                                                                                                     |                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                |

| Guidance and manufacturer's declaration - electromagnetic immunity                                                                                                                |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The AC-80MD is intended for use in the electromagnetic environment specified below. The customer or the user of the AC-80MD should assure that it is used in such as environment. |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Immunity test                                                                                                                                                                     | IEC 60601 test level        | Compliance level | Electromagnetic environment-guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Conducted RF<br>IEC 61000-4-6                                                                                                                                                     | 3 Vrms<br>150 kHz to 80 MHz | 3 Vrms           | <p>Portable and mobile RF communications equipment should be used no closer to any part of the AC-80MD, including cables, than the recommended separation distance calculated from the equation appliance to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1.2 \sqrt{P}$                                                                                                                                                                                                                                                                                                                                                                                              |
| Radiated RF<br>IEC 61000-4-3                                                                                                                                                      | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            | <p><math>d = 1.2 \sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.3 \sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<br/><sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with following symbol:</p>  |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AC-80MD is used exceeds the applicable RF compliance level above, the AC-80MD should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the AC-80MD.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

| Recommended separation distances between portable and mobile RF communications equipment and the AC-80MD                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                |                                         |                                          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------|------------------------------------------|
| The AC-80MD is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AC-80MD can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AC-80MD as recommended below, according to the maximum output power of the communications equipment.                                                                                                                                                                                                                 |                                                                |                                         |                                          |
| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Separation distance according to frequency of transmitter<br>m |                                         |                                          |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | 150 kHz to 80 MHz<br>$d = 1.2 \sqrt{P}$                        | 80 MHz to 800 MHz<br>$d = 1.2 \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.3 \sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 0.12                                                           | 0.12                                    | 0.23                                     |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 0.38                                                           | 0.38                                    | 0.73                                     |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 1.2                                                            | 1.2                                     | 2.3                                      |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | 3.8                                                            | 3.8                                     | 7.3                                      |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 12                                                             | 12                                      | 23                                       |
| <p>For transmitters rated a maximum output power not listed above, the recommended separation distance <math>d</math> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> |                                                                |                                         |                                          |

## Caution

When you dispose of the unit or accessories, you must obey the laws in the relative area or country and the regulations in the relative hospital regarding environmental pollution.

## Caution

The apparatus shall not be exposed to dripping or splashing. No objects filled with liquids, such as vases, shall be placed on the apparatus.

## Warning on power connection

Use a proper power cord for your local power supply.

1. Use the approved Power Cord (3-core mains lead) / Appliance Connector / Plug with earthing-contacts that conforms to the safety regulations of each country if applicable.
2. Use the Power Cord (3-core mains lead) / Appliance Connector / Plug conforming to the proper ratings (Voltage, Ampere).

If you have questions on the use of the above Power Cord / Appliance Connector / Plug, please consult a qualified service personnel.

## Warning on power connection for medical use

Please use the following power supply cord. With connectors (plug or female) and cord types other than those indicated in this table, use the power supply cord that is approved for use in your area.

|                                                | United States and Canada   |
|------------------------------------------------|----------------------------|
| Plug Type                                      | HOSPITAL GRADE*            |
| Cord Type                                      | Min.Type SJT<br>Min.18 AWG |
| Minimum Rating for Plug and Appliance Couplers | 10A/125V                   |
| Safety Approval                                | UL Listed and CSA          |

\*Note: Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked 'Hospital Only' or 'Hospital Grade'.

## For the customers in Europe

This product has been manufactured by or on behalf of Sony Corporation, 1-7-1 Konan Minato-ku Tokyo, 108-0075 Japan.

Inquiries related to product compliance based on European Union legislation shall be addressed to the authorized representative, Sony Deutschland GmbH, Hedelfinger Strasse 61, 70327 Stuttgart, Germany.

For any service or guarantee matters, please refer to the addresses provided in the separate service or guarantee documents.

## For the customers in the U.S.A.

**SONY LIMITED WARRANTY** - Please visit <http://www.sony.com/psa/warranty> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

## For the customers in Canada

**SONY LIMITED WARRANTY** - Please visit <http://www.sonybiz.ca/solutions/Support.do> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

## For the customers in Europe

Sony Professional Solutions Europe - Standard Warranty and Exceptions on Standard Warranty.

Please visit <http://www.pro.sony.eu/warranty> for important information and complete terms and conditions.

## For the customers in Korea

**SONY LIMITED WARRANTY** - Please visit <http://bpeng.sony.co.kr/handler/BPAS-Start> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

## Overview


This AC adaptor is designed for use with Sony medical equipment only.

*For details on connecting and using this unit, refer to the manual for your equipment.*

### Note

Do not install this equipment in a confined space, such as a bookshelf or similar unit.

## Location and Function of Parts and Controls

### 1 LED lamp

Lights when the power is on.

### 2 AC IN connector

Connect the AC power cord.

### 3 DC connector

Connect the DC cable to the DC connector of the medical equipment.

### Caution

Connect the DC connector to the equipment, and then connect the AC power cord of the AC adaptor.

When disconnecting the DC connector, disconnect the AC power cord of the AC adaptor before disconnecting the DC connector.


Do not touch the patient and the pin of the DC connector at the same time.

The pin of the DC connector applies a voltage of 24 V, which may harm the patient.

---

# How to Connect the DC Cable

Insert the DC cable into the DC connector of the medical equipment with its arrow mark facing upwards.


Make sure that the DC cable is connected securely until it is locked.

---

## Disconnecting the DC Cable

Hold the DC connector to pull the DC cable out.

Pull the connector towards you to release the lock.

### Note

Before connecting or disconnecting the DC cable, be sure to unplug the power plug from the mains outlet.

---

# Maintenance

---

## Cleaning the AC Adaptor

Use a wrung out damp cloth, or a damp cloth soaked in water and a mild detergent, to clean the AC adaptor surface.

To avoid damaging the AC adaptor, never use alcohol, thinner, or any other type of organic solvent, or any type of abrasive cleaner.

# Specifications

|                                   |                                                                                                                                                                                                                                                                                                                                                                                                       |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power requirement                 | AC 100 V - 240 V, 50/60 Hz                                                                                                                                                                                                                                                                                                                                                                            |
| Input current                     | 1.0 A - 0.5 A                                                                                                                                                                                                                                                                                                                                                                                         |
| Output voltage                    | DC 24 V                                                                                                                                                                                                                                                                                                                                                                                               |
| Output current                    | 3.3 A Max.                                                                                                                                                                                                                                                                                                                                                                                            |
| Operating temperature             | 0 °C to 40 °C (32 °F to 104 °F)                                                                                                                                                                                                                                                                                                                                                                       |
| Operating humidity                | 20% to 80%                                                                                                                                                                                                                                                                                                                                                                                            |
| Operating pressure                | 700 hPa to 1,060 hPa                                                                                                                                                                                                                                                                                                                                                                                  |
| Storage and transport temperature | -20 °C to +60 °C (-4 °F to +140 °F)                                                                                                                                                                                                                                                                                                                                                                   |
| Storage and transport humidity    | 20% to 90%                                                                                                                                                                                                                                                                                                                                                                                            |
| Storage and transport pressure    | 700 hPa to 1,060 hPa                                                                                                                                                                                                                                                                                                                                                                                  |
| Dimensions                        | 63 × 30.5 × 140 mm (w/h/d)<br>(2 1/2 × 1 1/4 × 5 5/8 inches) Not<br>incl. projecting parts                                                                                                                                                                                                                                                                                                            |
| Mass                              | 350 g (12 oz) (AC adaptor only)                                                                                                                                                                                                                                                                                                                                                                       |
| Accessories provided              | Instructions for Use (1)<br>Service Contact List (1)                                                                                                                                                                                                                                                                                                                                                  |
| Medical Specifications            | Protection against electric shock:<br>Class I<br>Protection against harmful ingress<br>of water:<br>Ordinary<br>Degree of safety in the presence of<br>a flammable anesthetic mixture<br>with air or with oxygen or<br>nitrous oxide:<br>Not suitable for use in the<br>presence of a flammable<br>anesthetic mixture with air or<br>with oxygen or nitrous oxide<br>Mode of operation:<br>Continuous |

The design and these specifications are subject to change without prior notification.

## Notes

- Always verify that the unit is operating properly before use. SONY WILL NOT BE LIABLE FOR DAMAGES OF ANY KIND INCLUDING, BUT NOT LIMITED TO, COMPENSATION OR REIMBURSEMENT ON ACCOUNT OF THE LOSS OF PRESENT OR PROSPECTIVE PROFITS DUE TO FAILURE OF THIS UNIT, EITHER DURING THE WARRANTY PERIOD OR AFTER EXPIRATION OF THE WARRANTY, OR FOR ANY OTHER REASON WHATSOEVER.
- SONY WILL NOT BE LIABLE FOR CLAIMS OF ANY KIND MADE BY USERS OF THIS UNIT OR MADE BY THIRD PARTIES.
- SONY WILL NOT BE LIABLE FOR THE TERMINATION OR DISCONTINUATION OF ANY SERVICES RELATED TO THIS UNIT THAT MAY RESULT DUE TO CIRCUMSTANCES OF ANY KIND.

If the unit is suddenly taken from a cold to a warm location, or if ambient temperature suddenly rises, moisture may form on the outer surface of the unit and/or inside of the unit. This is known as condensation. If condensation occurs, turn off the unit and wait until the condensation clears before operating the unit. Operating the unit while condensation is present may damage the unit.

The life expectancy of the AC adapter and the electrolytic capacitor is about 5 years under normal operating temperatures and normal usage (8 hours per day; 25 days per month). If usage exceeds the above normal usage frequency, the life expectancy may be reduced correspondingly.

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© 2013 Sony Corporation